

## Setting up Diasend to connect to Dexcom G5

1. Go to [www.diasend.com](http://www.diasend.com)
2. Create account in child's name
3. When prompted, enter **hospital ID 45-33954** (this will allow us to see the information)
4. Download the appropriate downloader software - Windows or Mac (see screenshot below)
5. Go to 'Account admin', then 'Connect app' tab (see screenshot below)
6. Click on 'Dexcom G5 Mobile'
7. Enter Dexcon login details.
8. You will start to see Dexcom blood sugars within five hours.

Glooko and Diasend have merged to make diabetes management easier. For more information about the diasend® product, please visit [glooko.com/diasend](http://glooko.com/diasend)

Existing user? Please login as usual.

Want to register a new patient account? Please click the yellow “Register here” button.

### Recent updates

July 2017.pdf

### Help Center

In need of assistance, tutorials or documentation? Click on one of the icons below.



UK - English 

**Username**

Enter Username

**Password**

Enter Password

Login

Forgot your password?

[Register here](#)

### Ready for a new career?

Click here to check out our open positions.

Sweden

USA

1

Login details

2

Personal information

3

Share data

4

Confirm registration

## Thank you

Your diasend® Personal account has now been registered.

In order to upload your data to your diasend® Personal account, you will need the tool diasend® Uploader (free of charge) installed on your computer.



After installing the diasend® Uploader application on your PC or Mac, you will be able to upload your data to your diasend® account. To download instructions regarding the installation and use of diasend® Uploader, please go to the support tab on our webpage. Please note: You may require a downloading cable specific to your device in order to upload your data. Please refer to the user guide of your device or contact your device vendor.

Thank you for choosing diasend®! You may now [log in!](#)

[Share data](#)[Connect app](#)[Recent Updates \(May 2017\)](#)[Customer Support](#)

## John Test Test

PDF

[Share data](#)[Export data](#)[Subscription](#)[Browse & connect apps](#)

IMPORTANT: It may take a couple of hours or until next daybreak before data may appear in your diasend® reports after you have established connection. Depending on activity tracker, your activity tracker's data is then automatically pushed to diasend® at least once a day when new data is available.



Fitbit

Fitbit offers wearable devices to help you lead a healthier, more active life.

[Connect](#)

Jawbone UP

UP is a wristband and app that tracks how you sleep, move and eat.

[Connect](#)

Moves App

Moves App is an activity diary of your life.

[Connect](#)

Runkeeper

The Runkeeper app is a personal trainer in your pocket.

[Connect](#)

Dexcom G5 Mobile

The Dexcom G5 Mobile is a sensor and app that continuously tracks your glucose levels and trends for better diabetes management.

[Connect](#)

Abbott FreeStyle LibreLink

The Freestyle Libre sensor is read via the LibreLink app in your Android smart phone. Connecting your LibreLink account to your diasend account means that your data will continuously sync with diasend.

[Connect](#)

Eversense

The Eversense Mobile Application works with the Eversense Continuous Glucose Monitoring System. The mobile app provides easy access to CGM data and eliminates the need to carry a dedicated receiver device.

[Connect](#)

# Your Dexcom Account

Sign in to manage your account.



Username \*

Password \*

Sign in >

[Forgot your Dexcom username or password?](#)

# Need a Dexcom Account?

You need a Dexcom Account to access mobile apps and web applications.

Create a Dexcom Account >