

Setting up Diasend to connect to Libre

1. Go to www.diasend.com
2. Create account in child's name
3. When prompted, enter **hospital ID 45-33954** (this will allow us to see the information)
4. Download the appropriate downloader software - Windows or Mac (see screenshot below)
5. Go to 'Account admin', then 'Connect app' tab (see screenshot below)
6. Click the LibreLink option
7. You should then be asked to connect the Libre with the yellow cable.

Glooko and Diasend have merged to make diabetes management easier. For more information about the diasend® product, please visit glooko.com/diasend

Existing user? Please login as usual.

Want to register a new patient account? Please click the yellow “Register here” button.

Recent updates

July 2017.pdf

Help Center

In need of assistance, tutorials or documentation? Click on one of the icons below.



UK - English

Username

Enter Username

Password

Enter Password

Login

Forgot your password?

Register here

Ready for a new career?

Click here to check out our open positions.

Sweden

USA

1

Login details

2

Personal information

3

Share data


4

Confirm registration

Thank you

Your diasend® Personal account has now been registered.

In order to upload your data to your diasend® Personal account, you will need the tool diasend® Uploader (free of charge) installed on your computer.



After installing the diasend® Uploader application on your PC or Mac, you will be able to upload your data to your diasend® account. To download instructions regarding the installation and use of diasend® Uploader, please go to the support tab on our webpage. Please note: You may require a downloading cable specific to your device in order to upload your data. Please refer to the user guide of your device or contact your device vendor.

Thank you for choosing diasend®! You may now [log in!](#)

Connect app

Recent Updates (May 2017)

Customer Support

John Test Test

PDF

Share data

Export data

Subscription

Browse & connect apps

IMPORTANT: It may take a couple of hours or until next daybreak before data may appear in your diasend® reports after you have established connection. Depending on activity tracker, your activity tracker's data is then automatically pushed to diasend® at least once a day when new data is available.



Fitbit

Fitbit offers wearable devices to help you lead a healthier, more active life.

[Connect](#)



Moves App

Moves App is an activity diary of your life.

[Connect](#)



Dexcom G5 Mobile

The Dexcom G5 Mobile is a sensor and app that continuously tracks your glucose levels and trends for better diabetes management.

[Connect](#)



Eversense

The Eversense Mobile Application works with the Eversense Continuous Glucose Monitoring System. The mobile app provides easy access to CGM data and eliminates the need to carry a dedicated receiver device.

[Connect](#)



Jawbone UP

UP is a wristband and app that tracks how you sleep, move and eat.

[Connect](#)



Runkeeper

The Runkeeper app is a personal trainer in your pocket.

[Connect](#)



Abbott FreeStyle LibreLink

The Freestyle Libre sensor is read via the LibreLink app in your Android smart phone. Connecting your LibreLink account to your diasend account means that your data will continuously sync with diasend.

[Connect](#)