

Somnotouch safety checklist

To be completed by parent or guardian on receipt of equipment

This is a safety checklist for the Somnotouch which has been couriered to you for use with your child. When you are ready to start the overnight sleep study, please follow the user guide for this piece of equipment.

1. Is the equipment dry? There should be no moisture in the bag surrounding the equipment.
2. Is the blue effort band adapter cable intact? Check that there is no damage to the outer cable sheath. There should be no internal cables or wires visible.
3. Check that the outer casing of the Somnotouch is intact and free from cracks or breakages.

If your child has been sent a camera in the box of equipment, please complete the following checks:

1. Check that the camera does not rattle when picked up and moved around.
2. Check that the camera power cable is intact and that there are no internal cable or wires visible.
3. Check that the camera power cable plug is free from cracks and is intact.

If you notice that anything is damaged please contact the sleep team immediately on 023 8120 8363 before any equipment is used.

For a translation of this document, or a version in another format such as easy read, large print, Braille or audio, please telephone **023 8120 4688**.

For help preparing for your visit, arranging an interpreter or accessing the hospital please visit **www.uhs.nhs.uk/additionalneeds**

