

Videosomnography safety checklist

To be completed by parent or guardian on receipt of equipment

This is a safety checklist for the videosomnography which has been couriered to you for use with your child. When you are ready to start the overnight sleep study, please follow the user guide for this piece of equipment.

Equipment

- Check that the equipment is dry. There should be no moisture in the bag surrounding the equipment.
- Check that the watch is intact. There should be no cables exposed, and the watch LED graphics should be clear and undamaged.

Video camera

- Check that the camera does not rattle when picked up and moved around.
- Check that the camera power cable is intact and that there are no internal cables or wires visible.
- Check that the camera power cable plug is intact and free from cracks.

Tripod

- Check that the tripod is stable and that there are no breakages.
- Check that the clip on top of the tripod is working and the camera fits into it.

Contact us

If you notice that anything is damaged, please contact us immediately before using any of the equipment.

Sleep team

Telephone: **023 8120 8363**

For a translation of this document, or a version in another format such as easy read, large print, Braille or audio, please telephone **0800 484 0135** or email **patientsupporthub@uhs.nhs.uk**

For help preparing for your visit, arranging an interpreter or accessing the hospital please visit **www.uhs.nhs.uk/additionalsupport**

www.uhs.nhs.uk/childrenshospital

Version 1. Published March 2021. Due for review March 2024. 2887